



Cold Fusion Cloning Kit

Cat. #s MC010A-1, MC100A-1, MC101A-1

User Manual

**Store the master mixture and positive controls at -20°C
Store the competent cells at -80°C.**

A limited-use label license covers this product. By use of this product, you accept the terms and conditions outlined in the Licensing and Warranty Statement contained in this user manual.

Contents

I.	Introduction	2
A.	Key Features	2
B.	Applications	2
C.	List of Components	3
D.	Storage	4
E.	Other Reagents Needed	4
F.	Technical Information	4
II.	Protocol	5
A.	Overview	5
B.	Preparation of Linearized Vector	5
C.	Primer Design	6
D.	PCR Requirements	6
E.	Cold Fusion Reaction	6
F.	Transformation	7
III.	Examples	8
A.	Primer Design for Positive Control Reaction	8
B.	Cloning a single DNA Fragment	8
C.	Joining Multiple DNA Fragments	8
IV.	Troubleshooting	9
V.	Technical Support	11
VI.	Licensing and Warranty	11

I. Introduction

The Cold Fusion technology is a simple, rapid and highly efficient PCR cloning kit. It allows you to directly clone any PCR product(s) to any linearized expression vector, at any site. The PCR fragments can be generated by Taq DNA polymerase or other high fidelity DNA polymerases, with primers that are designed to have ~15 bases of homology at the linear ends where the DNA-of-interest will “fuse”. The linearized vector can be generated by PCR or restriction enzyme digest (single or double cut). With a one tube simple reaction format, a 5 minute incubation at room temperature, followed by 10 minutes on ice, your PCR product(s) rapidly and accurately fuse into the linearized vector in the desired orientation. The kit is so robust that multiple DNA fragments can be assembled simultaneously and cloned into one construct in a single step. The system is highly efficient, with more than a 95% positive cloning rate.

A. Key Features

- Cloning is simple, rapid, accurate and directional
- Clone any insert, at any site within any vector
- Restriction enzyme, phosphatase and ligase free system
- Broad PCR size
- Joining multiple fragments at once
- High efficiency with > 95% positive clones

B. Applications

- PCR cloning into any vector
- Gene transfer from one vector to another
- Add adaptor, linker and tag before or after the insert
- Gene synthesis
- High throughput cloning

Cold Fusion™ Cloning Kit
Cat. # MC010A-1, MC 100A-1, MC 101A-1

C. List of Components

Cat. No. MC010A-1

5x Master Mixture (20rxns)	20µl
Linearized vector, positive control	5µl
500bp PCR insert, positive control	5µl
Competent cells (1x10 ⁹ cfu/µg)	10 tubes
User manual	1

Cat. No. MC100A-1

5x Master Mixture (20rxns)	40µl
Linearized vector, positive control	5µl
500bp PCR insert, positive control	5µl
Competent cells (1x10 ⁹ cfu/µg)	20 tubes
User manual	1

Cat. No. MC101A-1

5x Master Mixture (50rxns)	100µl
Linearized vector, positive control	10µl
500bp PCR insert, positive control	10µl
Competent cells (1x10 ⁹ cfu/µg)	50 tubes
User manual	1

D. Storage

Store the master mixture and positive controls at -20°C. Store the competent cells at -80°C.

E. Other Reagents Needed

- Gene-specific primers
- dNTPs, Taq or other high fidelity polymerase, and corresponding buffers for PCR
- QIAquick PCR Purification Kit (Cat # 28106, Qiagen)
- QIAquick Gel Extraction kit (Cat # 28704, Qiagen)
- SOC or LB Broth for transformation of bacteria
- LB + 100 µg/ ml Ampicillin plates

F. Technical Information


The competent cells provided with the Cold Fusion Cloning Kit has the following genotype:

F' {proAB+ lacIq lacZ Δ M15 Tn10(TetR) Δ (ccdAB)} mcrA Δ (mrr-hsdRMS-mcrBC) ϕ 80(lacZ) Δ M15 Δ (lacZYA-argF) U169 endA1 recA1 supE44 thi-1 gyrA96 relA1 tonA panD

Due to the presence of the Tn10 transposon element, which encodes the tetracycline-resistance gene (TetR), the competent cell will be able to grow in LB-Agar plates containing 10 µg/ml of tetracycline. Cloning of inserts into plasmids containing a TetR marker and transforming the competent cells in the kit will likely lead to the failure of the cloning reaction. In this case, we would recommend another high-quality competent cell that does not have built-in tetracycline resistance (TetR) for best results.

II. Protocol

A. Overview



B. Preparation of Linearized Vector

Complete linearization of the vector is critical to achieve a successful Cold Fusion cloning reaction. Incomplete linearization of the vector will result in high background. The linearized vector can be generated by PCR or restriction enzyme digest (single or double digest) and should be purified using a gel or PCR purification kit.

Due to the digestion efficiency, different restriction enzymes will generate different levels of background. In general, two enzyme digestion is better than a single enzyme digestion. The further the restriction sites are apart, the better the digestion efficiency. Increasing the enzyme digestion time and the digestion reaction volume will also help reduce the background. For many enzymes, we recommend incubate the digestion reaction between 3 hours and overnight in order to increase linearization and reduce background.


Check the background of your vector by transforming 1 μ g (10-100ng) linearized and purified vector into competent cells. If the background is high, continue digesting the remaining vector for a longer time after the addition of more restriction enzyme(s).

We recommend digesting 2 μ g vector in 50 μ l reaction overnight. Use QIAGEN's QIAquick Spin Gel Extraction kit for gel purification and elute the DNA with 30 μ l dH₂O.

C. Primer Design


Forward Primer Design

The forward primer should contain 18-20 bases complimentary to the 5' end of your gene-of-interest plus 15 bases corresponding to the vector. Depending on the purpose of your cloning, your primers might also need to contain a Kozak sequence and/ or ATG start site.



Reverse Primer Design

Reverse primers should be made to the negative strand of DNA. Depending on the vector you are using, the 3' end of the insert may need to contain a stop codon (TAG, TAA, or TGA). The 3' primer should contain 18-20 bases complimentary to the 3' end of your gene-of-insert, plus 15 bases corresponding to the vector.



Introducing Restriction Sites

A restriction site can be introduced in the middle of the primer and can be the same as or different than the one used to linearize the vector. For multiple DNA fragment joining, it is recommended that each PCR product shares ~ 18 base pairs of homology.

D. PCR Requirements

- The PCR fragments for the cDNA-of-interest should be generated using *Taq* DNA polymerase or other high fidelity DNA polymerase.
- The melting temperature (T_m) should be calculated based on the gene-specific ends of the primer, NOT the entire primer.
- Specific PCR reaction conditions should be optimized for the cDNA-of-interest.
- After completion of the PCR reaction, gel purify the appropriate band to remove any extra primers or primer dimers that will inhibit the Cold Fusion reaction. We recommend the QIAquick PCR Purification Kit (Cat # 28106, Qiagen).

E. Cold Fusion Reaction

Set up the following reaction in a 1.5 ml sterile reaction tube by mixing the following reagents gently and then spin down briefly to collect the reagents at the bottom of the tube.

Cloning reaction

Linearized destination vector (10-100ng/ μ l)	1 μ l*
PCR insert(s) (20-200ng/ μ l) (for each PCR Product)	1 μ l*
dH ₂ O	_ μ l

Cold Fusion™ Cloning Kit
Cat. # MC010A-1, MC 100A-1, MC 101A-1

5x master mix	2µl
total	10µl

Positive control reaction

Linearized vector (positive control)	1µl
500bp PCR insert (positive control)	1µl
dH ₂ O	6µl
5x master mix	2µl
total	10µl

Negative Control

Linearized destination vector (10-100ng/µl)	1µl*
dH ₂ O	7µl
5x master mix	2µl
total	10µl

* 2:1 or 1:1 molar ratio of insert: vector works well in the Cold Fusion reaction.

For reactions with larger volumes of vector and insert (>8µl of vector + insert), double the amount of reaction buffer and enzyme, and add dH₂O for a total volume of 20µl.

When using Cold Fusion cloning kit for the first time, we strongly recommend that you perform the positive and negative control reactions in parallel with your Cold Fusion cloning reaction. The positive control 500bp PCR insert and linearized vector provided in the kit have already been purified. There is no treatment needed prior to the cloning reaction.

Cold Fusion Reaction Incubation

1. 5 minutes at room temperature
2. 10 minutes on ice

F. Transformation

1. Add 50µl Cold Fusion competent cells to the cloning mixture
2. Incubate on ice for 20 minutes
3. Heat shock at 42°C for 50 seconds
4. Transfer on ice for 2 minutes
5. Add 250µl S.O.C medium or LB broth
6. Incubate at 37°C for an hour
7. Take 100µl culture spread on pre-warmed (37°C) culture plate containing Ampicillin
8. Incubate the plate at 37°C overnight


III. Examples

A. Primer Design for Positive Control Reaction

1 Original Vector Sequence. Linearize with EcoRI and BamHI

5'- ACGACGTTGAAAACGACGGCGAGTGAATTCGAGCTCGGTACCCGGGATCCTCTAGAGTCGACCTGCAGGCATGCAAGCT -3'

2 Construct-of-interest shown below. Black = Vector sequence. Blue = Insert sequence. Red = restriction site




3 Forward primer Sequence: 5'-AAAACGACGGCCAGTGAATTCatggagagcgcacgagcggc-3'


Reverse primer Sequence: 5'- CAGTCGACTCTAGAGGATCCcttgaagtgcattgtggctgc -3'

B. Cloning a single DNA Fragment

Broad PCR product cloning range with high efficiency




Good for any vector with 5' overhangs, 3' overhangs and blunt ends



C. Joining Multiple DNA Fragments

Note: For multiple DNA fragments cloning, depending on the number and the size of each insert, you may obtain fewer colonies than those with one or two fragments.

Assemble and clone multiple DNA fragments simultaneously into one construct in a single step



IV. Troubleshooting

Problem	Probable cause	Solution
1. No or few colonies obtained from the transformation	Primer sequences are incorrect	Check primer sequences to ensure that they provide 15 bases of homology with the region flanking the insertion site.
	Suboptimal PCR product	Optimize your PCR amplification reactions so that you generate pure PCR products. Use a different method to purify your PCR product.
	Low DNA concentration in reaction	It is imperative to obtain as high a DNA concentration as possible in your Cold Fusion reaction.
	Inhibitory contaminants from PCR product or linearized vector	Both the PCR product and the linearized vector should be purified.
	Transform with too much reaction mixture	Do not add more than 10µl of reaction mixture to 50µl of competent cells. Too much reaction mixture inhibits the transformation.
	Low quality or poor handling of competent cells	Handle the competent cells gently. Do not re-freeze cells after thawing. Quality of competent cells may be tested by transforming a circular plasmid to determine cells' competency. Competent cells with a transformation efficiency of 1×10^9 cfu/ µg are recommended.
	Wrong antibiotic or too much	Choose plates with the appropriate concentration of

Problem	Probable cause	Solution
	antibiotic in the media	the right antibiotic.
2. Large numbers of colonies contain no insert	Incomplete linearization of your vector	It is critical to remove any uncut vector prior to use in the Cold Fusion reaction. If necessary, re-digest your vector and gel purify.
	Contamination of cloning reaction with plasmid with the same antibiotic resistance	If you insert was amplified from a plasmid, circular DNA may have carried through purification and contaminated the cloning reaction. We recommended gel purifying your PCR product or linearizing the template DNA before performing PCR.
2. Large numbers of colonies contain no insert	Plates are too old or contained incorrect antibiotic	Make sure that your antibiotic plates are fresh. Check the antibiotic resistance of your fragment.
3. Clones contain incorrect insert	PCR products contain non-specifically amplified artifacts	Optimize your PCR reaction to improve the specificity. Screen more colonies for the correct clones.

V. Technical Support

For more information about SBI products and to download manuals in PDF format, please visit our web site:

<http://www.systembio.com>

For additional information or technical assistance, please call or email us at:

Phone: (650) 968-2200

(888) 266-5066 (Toll Free)

Fax: (650) 968-2277

E-mail:

General Information: info@systembio.com

Technical Support: tech@systembio.com

Ordering Information: orders@systembio.com

System Biosciences (SBI)
265 North Whisman Rd.
Mountain View, CA 94043

VI. Licensing and Warranty

Limited Use License

Use of the Cold Fusion™ Kit (i.e., the “Product”) is subject to the following terms and conditions. If the terms and conditions are not acceptable, return all components of the Product to System Biosciences (SBI) within 7 calendar days. Purchase and use of any part of the Product constitutes acceptance of the above terms. Purchase of the product does not grant any rights or license for use other than those explicitly listed in this Licensing and Warranty Statement. Use of the Product for any use other than described expressly herein may be covered by patents or subject to rights other than those mentioned. SBI disclaims any and all responsibility for injury or damage which may be caused by the failure of the buyer or any other person to use the Product in accordance with the terms and conditions outlined herein. SBI has pending patent applications related to the Product. For information concerning licenses for commercial use, contact SBI.

Limited Warranty

SBI warrants that the Product meets the specifications described in the accompanying Product Analysis Certificate. If it is proven to the satisfaction of SBI that the Product fails to meet these specifications, SBI will replace the Product or provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to SBI within 30 days of receipt of the Product.

SBI’s liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. SBI’s liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. SBI does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

SBI is committed to providing our customers with high-quality products. If you should have any questions or concerns about any SBI products, please contact us at (888) 266-5066.

© 2013 System Biosciences (SBI), All Rights Reserved.