

phiC31 Integrase Vector System


One-step gene addition technology



Nonviral transgene delivery

The phiC31 integrase is a sequence-specific recombinase encoded within the genome of the bacteriophage phiC31. The phiC31 integrase mediates recombination between two 34 base pair sequences termed attachment sites (*att*), one found in the phage and the other in the bacterial host. This serine integrase has been shown to function efficiently in many different cell types, including mammalian cells. In the presence of phiC31 integrase, an *attB*-containing donor plasmid can be unidirectionally integrated into a target genome through recombination at sites with sequence similarity to the native *attP* site (termed pseudo-*attP* sites). phiC31 integrase can integrate a plasmid of any size, as a single copy, and requires no cofactors. The integrated transgenes are stably expressed and heritable.


How phiC31 integrase works



Highlights

- Nonviral transgenesis
- Single-copy integration
- Preferred integration at active transcription sites
- Unlimited cargo capacity
- Easily make stable cell lines using transfection-based methods
- Suitable for engineering Stem Cells and In vivo animal models

phiC31 integrase donor vector formats




phiC31 Integrase Vector System

One-step phiC31-mediated gene addition

HEK293 cells seeded onto 6-well plates were transfected the next day in triplicate according to the manufacturer's instructions with Lipofectamine 2000 (Invitrogen) and 40 ng pFC-PGK-MCS-EF1-GFP-Puro Dual Promoter phiC31 Donor Vector (Cat# FC551A-1) and either 2 µg carrier plasmid DNA or 2 µg phiC31 integrase expression plasmid. The cells were trypsinized after 48 hours, resuspended in 1 ml of media, and split 1:20 onto 10 cm plates containing 10 ml of media. After 24 hours, puromycin 0.5 µg/ml was added to the medium to begin selection. After 11 days, the cells were fixed and stained with a solution of 50% methanol plus 1% methylene blue. The plates were washed twice with PBS and allowed to air dry. The number of visible colonies were imaged and colony number counts assessed. There were no visible colonies on the minus phiC31 integrase and over 300 colonies were counted on the plus phiC31 integrase plate (colony assay data shown to the right).


- phiC31
Integrase

+ phiC31
Integrase




phiC31 donor integration with GFP and Puro marker expression

Approximately 400,000 HEK293 cells were co-transfected with 1.9 µg Integrase expression vector with 40ng pFC-PGK-MCS-EF1-GFP-Puro Dual Promoter phiC31 Donor Vector (Cat# FC551A-1). The cells were split 1:10 after 24 hours post-transfection. Puromycin selection at 1 µg/ml was initiated after another 24 hours with continuous selection for 4 days. Cells were imaged at 4 and 11 days. Representative GFP fluorescence and phase contrast fields are shown below.



4 Days

11 Days



We Also Offer Custom Services - Custom phiC31 donor vector cloning and cell lines

System Biosciences offers a wide-range of custom services to support your research, allowing you to spend less time making tools, and more time making discoveries. To learn more, visit our website at www.systembio.com/service or call us at 888-266-5066.



System Biosciences, Inc.
265 North Whisman Rd.
Mountain View, CA 94043

Toll Free: 888.266.5066
Fax: 650-968-2277
Email: info@systembio.com
www.systembio.com

© 2013 System Biosciences, Inc. All rights reserved. System Biosciences and the System Biosciences logo are trademarks of System Biosciences, Inc.