



SBI offers Express Lentiviral Packaging service to produce high-quality, high-titer virus particles containing your construct in only 2-3 weeks turnaround time.

Choose your level of titer:

Packaging Service	Catalog #	Titer	IFUs	Price
Regular Titer (RT)	CS960A-1	>10 <sup>7</sup> IFU/mL	2.5 x 10 <sup>6</sup>	\$1,360
High Titer (HT)	CS965A-1	>10 <sup>8</sup> IFU/mL	25 x 10 <sup>6</sup>	\$2,048
Ultra-High Titer (UHT)	CS970A-1	>10 <sup>9</sup> IFU/mL	250 x 10 <sup>6</sup>	\$2,573

**Service Details:** SBI's packaging service is available for 2nd and 3rd generation HIV- or FIV-based lentivectors only. The lentivector must be capable of producing virus, with no internal poly A signals, toxic gene load or unusual 2° structure. For SBI to titer your virus, your lentivector must contain a WPRE element. See the Plasmid Info Sheet for further details and requirements.

**You** provide SBI with endotoxin-free plasmid DNA of your lentivector (RT: 20 µg, HT: 40 µg, UHT: 80 µg), shipped on dry ice. Plasmid Info Sheet required.

**SBI** provides ten 25 µL aliquots of VSV-G pseudotyped lentiviral particles at the specified level of titer, along with titer measurement based on infectious particles.

Shipping Address:	Billing Address:
Name:	Name:
Institution:	Institution:
Mailing Address:	Mailing Address:
City:	City:
State:	State:
Zip:	Zip:
Telephone:	Telephone:
Fax:	Fax:

Payment by Purchase Order	Payment by Credit Card
PO number:	Card #:
Expiration date:	Mastercard or Visa
	Expiration Date:
Please fax a copy of PO with this form	Name on Card: